

DEVICE SPECS

ENDOPATH® ETS Articulating Linear Cutters



ENDOPATH™ ETS Articulating Linear Cutters

Codes	ATS45
Associated Reloads	TR45W, 6R45B
Description	45mm ETS Articulating Linear Cutter
Firing	Single Stroke
Articulation	Manual
Shaft Length	34cm
Degrees of Articulation	70°
Staple Rows	6 (when using TR45W or 6R45B)
Max Firings	8
Size	45mm
Staple Form	B-Form
Qty./Box	3