

**DEVICE SPECS**

# PROXIMATE™ Reloadable Linear Stapler (TX)



**PROXIMATE™ Reloadable Linear Stapler (TX)  
with Vascular Reload**



**PROXIMATE™ Reloadable Linear Stapler (TX)  
with Blue Reload**



**PROXIMATE™ Reloadable Linear Stapler (TX)  
with Green Reload**

<b>Codes</b>	TX30V	TX30B, TX60B	TX30G, TX60G
<b>Associated Reloads</b>	XR30V	XR30B, XR60B	XR30G, XR60G
<b>Description</b>	30mm Vascular Linear Stapler	30mm & 60mm Linear Stapler	30mm & 60mm Linear Stapler
<b>Device Technology</b>	Auto Retaining Pin Advancement	Auto Retaining Pin Advancement	Auto Retaining Pin Advancement
<b>Max Firings</b>	8	8	8
<b>Staple Rows</b>	3	2	2
<b>Staple Line Length</b>	30mm	30mm, 60mm	30mm, 60mm
<b>Closed Staple Height</b>	1.0mm	1.5mm	2.0mm
<b>Open Staple Height</b>	2.4mm	3.6mm	4.7mm
<b>Staple Technology</b>	B-Form	B-Form	B-Form
<b>Tissue Thickness</b>	Vascular	Regular	Thick
<b>Qty./Box</b>	3	3	3